Whitepaper DOAC Dipstick

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Inhalt:

GRUNDLEGENDE PUBLIKATIONEN:

- I. Accuracy of a Rapid Diagnostic Test for the Presence of Direct Oral Factor Xa or Thrombin Inhibitors in Urine – A Multicenter Trial (Thromb Haemost, 2020)
- II. Detection of Direct Oral Anticoagulants in Patient Urine Samples by Prototype and Commercial Test Strips for DOACs – A Systematic Review and Meta-analysis (TH Open, 2021)
- III. Performance Characteristics of DOAC Dipstick in Determining Direct Oral Anticoagulants in Urine (Clin Appl Thromb Haemost, 2021)
- IV. Algorithm for Rapid Exclusion of Clinically Relevant Plasma Levels of Direct Oral Anticoagulants in Patients Using the DOAC Dipstick: An Expert Consensus Paper (Thromb Haemost, 2024)

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- V. DOAC Dipstick testing can reliably exclude the presence of clinically relevant DOAC concentrations in circulation (Thromb Haemost, 2022)
- VI. Performance of a Qualitative Point-of-Care Strip Test to Detect DOAC Exposure at the Emergency Department: A Cohort-Type Cross-Sectional Diagnostic Accuracy Study (Thromb Haemost, 2022)
- VII. Assessment of Direct Oral Anticoagulant Status Using the DOASENSE Dipstick in Thrombolysis Eligible Patients With Stroke: Proof-of-Concept Study (Stroke, 2023)
- VIII. Evaluation of DOAC Dipstick Test for Detecting Direct Oral Anticoagulants in Urine Compared with a Clinically Relevant Plasma Threshold Concentration (Clinical and Applied Thrombosis/Hemostasis, 2022)
- IX. Comparison of the DOAC Dipstick Test on Urine Samples With Chromogenic Substrate Methods on Plasma Samples in Outpatients Treated With Direct Oral Anticoagulants (Clin Appl Thromb Hemost, 2023)
- X. The DOAC (direct oral anticoagulant) Dipstick reliably excludes residual FXa inhibitors levels in the preoperative setting (Thromb Res, 2024)

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I. Accuracy of a Rapid Diagnostic Test for the Presence of Direct Oral Factor Xa or Thrombin Inhibitors in Urine — A Multicenter Trial

Job Harenberg, Jan Beyer-Westendorf, Mark Crowther, Jonathan Douxfils, Ismail Elalamy, Peter Verhamme, Rupert Bauersachs, Svetlana Hetjens, Christel Weiss

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Online-Link zur Publikation: <https://doi.org/10.1055/s-0039-1700545>

II. Detection of Direct Oral Anticoagulants in Patient Urine Samples by Prototype and Commercial Test Strips for DOACs – A Systematic Review and Meta-analysis

Andrea Martini, Job Harenberg, Rupert Bauersachs, Jan Beyer-Westendorf, Mark Crowther, Jonathan Douxfils, Ismail Elalamy, Christel Weiss, Svetlana Hetjens

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Online-Link zur Publikation: <https://doi.org/10.1055/s-0041-1732437>

III. Performance Characteristics of DOAC Dipstick in Determining Direct Oral Anticoagulants in Urine

Job Harenberg, Andrea Martini, Shanshan Du, Sandra Krämer, Christel Weiss, Svetlana Hetjens

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IV. Algorithm for Rapid Exclusion of Clinically Relevant Plasma Levels of Direct Oral Anticoagulants in Patients Using the DOAC Dipstick: An Expert Consensus Paper

Harenberg J, Gosselin RC, Cuker A, Becattini, Pabinger I, Poli S, Weitz JI, Ageno W, Bauersachs R, Ćelap I, Choi P, Douketis J, Douxfils J, Elalamy I, Fareed J, Falanga A, Favalor E, Gerotziakas GT, Herkner H, Hetjens S, Klamroth R, Heubner L, Langer F, Lip GYH, Mac Grory B, Margetić S, Merrelaar A, Pikta M, Renne T, Schwameis M, Schulman S, Strbian D, Tafur A, Violi F, Vassart J, Walenga J, Weiss C

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
	Direct oral factor Xa inhibitors Mean (95% CI)	Direct oral thrombin inhibitor Mean (95% CI)
Sensitivity	97.8 (95.6 – 99.0)	98.3 (91.0 – 100)
NPV	86.6 (76.0 – 93.7)	99.6 (97.7 – 100)
PPV	87.2 (83.7 – 90.1)	73.4 (63.7 – 83.2)
Specificity	50.0 (40.2 – 59.0)	91.8 (87.7 – 94.9)

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V. DOAC Dipstick testing can reliably exclude the presence of clinically relevant DOAC concentrations in circulation

Sandra Margetić, Ivana Čelap, Arijana Lovrenčić Huzjan, Marijana Bosnar Puretić, Sandra Šupraha Goreta, Anesa Čajević Glojnarčić, Diana Delić Brkljačić, Pavao Mioč, Job Harenberg, Svetlana Hetjens, Christel Weiss

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VI. Performance of a Qualitative Point-of-Care Strip Test to Detect DOAC Exposure at the Emergency Department: A Cohort-Type Cross-Sectional Diagnostic Accuracy Study

Anne E. Merrelaar, Magdalena S. Bögl, Nina Buchtele, Marieke Merrelaar, Harald Herkner, Christian Schoergenhofer, Job Harenberg, Jonathan Douxfils, Romain Siriez, Bernd Jilma, Alexander O. Spiel, Michael Schwameis

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VII. Assessment of Direct Oral Anticoagulant Status Using the DOASENSE Dipstick in Thrombolysis Eligible Patients with Stroke: Proof-of-Concept Study

Peter Shuangyue Tan, Peter SW Park, Ross Cody, Tanya Frost, Bailey McNamara, Marija Borosak, Philip MC Choi
Stroke. 2023;54: e142-e144.

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VIII. Evaluation of DOAC Dipstick Test for Detecting Direct Oral Anticoagulants in Urine Compared with a Clinically Relevant Plasma Threshold Concentration

Lenna Örd, Toomas Marandi, Marit Märk, Leonid Raidjuk, Jelena Kostjuk, Valdas Banys, Karit Krause, Marika Pikta

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IX. Comparison of the DOAC Dipstick Test on Urine Samples With Chromogenic Substrate Methods on Plasma Samples in Outpatients Treated With Direct Oral Anticoagulants

Loula Papageorgiou, Svetlana Hetjens, Jawed Fareed, Sanny Auge, Laetitia Tredler, Job Harenberg, Christel Weiss, Ismail Elalamy, Grigorios T. Gerotziapas

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X. The DOAC (direct oral anticoagulant) Dipstick reliably excludes residual FXa inhibitors levels in the preoperative setting

Julie Vassart, Alfonso Tafur, Marisa Durante, Anne-Yseult Boucher, Laure Morimont, Jeanine M. Walenga, Job Harenberg, Jonathan Douxfils

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XII. DOAC Dipstick for detecting direct oral anticoagulants

National Institute for Health and Care Excellence (NICE) des UK NHS
Medtech innovation briefing [MIB248], Februar 2021

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Online-Link zur Publikation: <https://www.nice.org.uk/advice/mib248>

XIII. 2021 Australasian Anesthesia Blue Book

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Online-Link zur Publikation (siehe "Blue Book 2021"-Link auf dieser Webseite): [https://www.anzca.edu.au/safety-and-advocacy/advocacy/college-publications/australasian-anaesthesia-\(the-blue-book\)](https://www.anzca.edu.au/safety-and-advocacy/advocacy/college-publications/australasian-anaesthesia-(the-blue-book))

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Clinical guideline on reversal of direct oral anticoagulants in patients with life threatening bleeding

Grottke, Oliver; Afshari, Arash; Ahmed, Aamer; Arnaoutoglou, Eleni; Bolliger, Daniel; Fenger-Eriksen, Christian; von Heymann, Christian

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Online-Link zur Publikation: <https://doi.org/10.1097/EJA.0000000000001968>

SIEHE AUCH:

XVII: Clinical guideline on reversal of direct oral anticoagulants in patients with life threatening bleeding: use of a urine dipstick

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XVIII: Recent guideline from the British Society of Haematology on the measurement of direct oral anticoagulants lists the DOASENSE point-of-care test

Baker P, Platton S, Arachchillage DJ, Kitchen S, Patel J, Riat R, Gomez K.

Br J Haematol. 2024 Oct;205(4):1302-1318.

Online-Link zur Publikation: <https://onlinelibrary.wiley.com/doi/10.1111/bjh.19729>

UND

XIX: Urgent Reversal of Direct Oral Anticoagulants in Critical and Life-Threatening Bleeding: A Multidisciplinary Expert Consensus

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EU Certificate

for the assessment of the technical documentation

according to Regulation on In-vitro Diagnostic Medical Devices (EU) 2017/746 Annex IX Chapter II

As a Notified Body of the European Union, DEKRA Certification GmbH certifies, that the manufacturer **DOASENSE GmbH**

Single Registration Number (SRN): **DE-MF-000015367**
Waldhofer Straße 102, 68123 Heidelberg, Germany

that the technical documentation of the product(s) described in the annex complies with the provisions of the Regulation on In-Vitro Diagnostic Medical Devices (EU) 2017/746. The certificate is based on the results of the assessment of the technical documentation according to the Regulation on In-Vitro Diagnostic Medical Devices (EU) 2017/746 Annex IX Chapter II, which are recorded in the report referred to in the annex.

Product: **DOAC Dipstick**

EU Certificate no.: **51551-71-A0-00** Certificate valid from: **2025-02-21**
Certificate valid to: **2030-02-20**

Digitally signed by Andrea Schmidt
Date: 2025.02.21 14:00:43+0100

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Annex to the EU Certificate no. 51551-71-A0-00

This certificate covers the following device:

Class C, Near patient testing	
Basic UDI-DI: 426055266DipstickH3	Intended Purpose:
Device Name: DOAC Dipstick	The diagnostic test strip DOAC Dipstick is intended for qualitative detection of the absence or presence of direct oral anticoagulants (DOACs: Dabigatran, Apixaban, Edoxaban, and Rivaroxaban) in human urine by visual identification of colours. The DOAC Dipstick is an in vitro diagnostic test and can be used at the Point of Care (POCT / Near Patient Test) or in the laboratory. The DOAC Dipstick is intended for professional use only.
Model: 0001	
Type: W01010802 Clinical Chemistry, Rapid Test & Point of Care	
Conformity assessment procedure: Annex IX, Chapter II Section 4, and 5.1, and Chapter III	

Change to previous certificate: n.a.

Remark: For the placing on the market of the product(s) referred to above, an additional EU certificate for the assessment of the quality management system in accordance with Annex IX, Chapter I + III is required.

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EU Certificate

for the assessment of the quality management system

according to Regulation on In-vitro Diagnostic Medical Devices (EU) 2017/746 Annex IX Chapter I+III

As a Notified Body of the European Union, DEKRA Certification GmbH certifies, that the manufacturer **DOASENSE GmbH**

Single Registration Number (SRN): **DE-MF-000015367**
Waldhofer Straße 102, 68123 Heidelberg, Germany

applies a quality management system according to Annex IX Chapter I+III of the regulation on in-vitro diagnostic medical devices (EU) 2017/746 for the devices listed in the annex. This certificate is based on the assessments listed in No51551-00 and is only valid in conjunction with the successful completion of the annual surveillance audits.

EU Certificate no.: **51551-70-00-00** Certificate valid from: **2025-02-21**
Certificate valid to: **2028-04-04**

Digitally signed by Andrea Schmidt
Date: 2025.02.21 14:00:43+0100

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Annex to the EU Certificate no. 51551-70-00-00

Following device/device categories are included in this certificate:

Class C near-patient test devices
DOAC Dipstick

Intended Purpose:
The diagnostic test strip DOAC Dipstick is intended for qualitative detection of the absence or presence of direct oral anticoagulants (DOACs: Dabigatran, Apixaban, Edoxaban, and Rivaroxaban) in human urine by visual identification of colours. The DOAC Dipstick is an in vitro diagnostic test and can be used at the Point of Care (POCT / Near Patient Test) or in the laboratory. The DOAC Dipstick is intended for professional use only.

Basic-UDI-DI: 426055266DipstickH3

- Level 1: Class C for near-patient testing
- Level 2: IVR005 Devices intended to be used for monitoring of levels of medicinal products, substances or biological components
- Level 3: W01010802 Clinical Chemistry, Rapid Test + Point-of-Care-Test
- Level 4: 426055266DipstickH3
- Level 5: IVP 3002 In vitro diagnostic devices which require knowledge regarding biochemistry

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Walenga JM

Can We Improve on the Rapid Assessment of Clinically Relevant Levels of Direct Acting Oral Anticoagulants (DOAC)?
Clin Appl Thromb Hemost. 2022 Jan-Dec;28:10760296221096422. doi: 10.1177/10760296221096422.

Lippi G, Favaloro EJ

Pearls and Pitfalls in the Measurement of Direct Oral Anticoagulants. Semin Thromb Hemost. 2024 Nov;50(8):1114-1122.

Qiao J, Tran MH

Challenges to Laboratory Monitoring of Direct Oral Anticoagulants.

Clin Appl Thromb Hemost. 2024 Jan-Dec;30:10760296241241524. doi: 10.1177/10760296241241524.

Baker P, Platton S, Arachchillage DJ, Kitchen S, Patel J, Riat R, Gomez K; BSH Committee

Measurement of heparin, direct oral anti-coagulants and other non-coumarin anti-coagulants and their effects on haemostasis assays: A British Society for Haematology Guideline.

Br J Haematol. 2024 Oct;205(4):1302-1318.

De Simone B, Chouillard E, Podda M, Pararas N, de Carvalho Duarte G, Fugazzola P, Birindelli A, Coccolini F, Polistena A, Sibilla MG, Kruger V, Fraga GP, Montori G, Russo E, Pintar T, Ansaloni L, Avenia N, Di Saverio S, Leppäniemi A, Lauretta A, Sartelli M, Puzziello A, Carcoforo P, Agnoletti V, Bissoni L, Isik A, Kluger Y, Moore EE, Romeo OM, Abu-Zidan FM, Beka SG, Weber DG, Tan ECTH, Paolillo C, Cui Y, Kim F, Picetti E, Di Carlo I, Toro A, Sganga G, Sganga F, Testini M, Di Meo G, Kirkpatrick AW, Marzi I, deAngelis N, Kelly MD, Wani I, Sakakushev B, Bala M, Bonavina L, Galante JM, Shelat VG, Cobianchi L, Mas FD, Pikoulis M, Damaskos D, Coimbra R, Dhesis J, Hoffman MR, Stahel PF, Maier RV, Litvin A, Latifi R, Biffi WL, Catena F

The 2023 WSES guidelines on the management of trauma in elderly and frail patients.

World J Emerg Surg. 2024 May 31;19(1):18. doi: 10.1186/s13017-024-00537-8.

Schulman S, Bhagirath V, Chan N, Germini F, Ikesaka R, Iorio A, Mithoowani S, Weitz JI, Gross PL

Prothrombin complex concentrate for emergency surgery in patients on oral Xa-inhibitors.

J Thromb Haemost. 2024 Oct;22(10):2761-2766.

Weber C, Rigby A, Lip GYH

Thrombosis and Haemostasis 2024 Editors' Choice Papers.

Thromb Haemost. 2025 Jan;125(1):85-91.

Tripodi A, Marchetti M, Scalabrino E

Direct oral anticoagulants measurement: when is it necessary?

Blood Transfus. 2025 Jan;23(1):55-58.

Callaly E, Tjahyadi J, Yin V, Holland B, Dewey H, Borosak M, Choi PMC

Results of rapid plasma direct oral anticoagulant levels greatly expand stroke thrombolysis eligibility: evidence from real-world implementation.

Intern Med J. 2025 Apr;55(4):664-668.

Ten Cate H

The wonders of anticoagulation.

Front Cardiovasc Med. 2025 Jan 13;11:1517109. doi: 10.3389/fcvm.2024.1517109.

Menichelli D, Pannunzio A, Baldacci E, Cammisotto V, Castellani V, Mormile R, Palumbo IM, Chistolini A, Violi F, Harenberg J, Pastori D, Pignatelli P

Plasma Concentrations of Direct Oral Anticoagulants in Patients with Nonvalvular Atrial Fibrillation and Different Degrees of Obesity.

Clin Pharmacokinet. 2025 Mar;64(3):453-462.

Mithoowani S, Tan CW, Siegal DM

Measuring Direct Oral Anticoagulant (DOAC) Levels: Applications, Limitations, and Future Directions.

Int J Lab Hematol. 2025 Apr 22. doi: 10.1111/ijlh.14483.

Kumar, Mandeep, Oprea, Adriana D

Perioperative Anticoagulation: Considerations Impacting Outcomes.

International Anesthesiology Clinics 2025, 63(2):44-53.